Preface

Reminding Before your testing

Breath tester should be used only to give an indication of the possible presence of alcohol in the breath. You should not rely upon it as the sole basis to determine intoxication or whether it is safe to drive a vehicle, operate equipment, or engage in dangerous activities.

Everyone has different body responds to alcohol consumption and his testing result only as reference, not a subject standard for consequent decision.

The manufacture, importer or distributor takes no responsibility whatsoever for the use of this product for any reason. This product must not be used as a tool for determining whether a person is able to operate a motor vehicle or device legally or safely. The intake of any alcohol will impair reflexes and judgment to operate motor vehicle.

Before your testing, please view the manual carefully and follow the instruction strictly.

General Introduction

Breath alcohol tester is kind of tester designed to measure concentration result of the breathed alcohol in the human body. This device adopts advanced MP-3 alcohol sensor which has excellent sensitivity and reproducibility, fashion and portable design make it more convenient for personal use. When the alcohol content exceeds the preset limited level, this device would send an audio and video warning to remind your safety.

Main Feature:

- Advanced alcohol sensor
- Quick response
- SMD assembling, stable performance
- Smart MCU control
- Direct testing process LED indication
- Digital LED display with red words
- Portable and fashion design
- Audio warning beyond pre-set limit
- Sensor Fault Self checking
- Battery saved design, low voltage indication
- Portable bag (optional)
- Car power supply

Technical Data

- Sensor type: Semi-conductor alcoholic sensor
- Detection Range: 0,00 – 2.1‰ BAC
- Alarming Level: 0.1 mg/l
- Accuracy 0,01‰ BAC
- Response Time: <5s
- Warm-up Time: <20s
- Resume Time: <20s
- Working Voltage: DC3.0V (2 x AAA) Batteries;
- Working Current: <140mA
- Best Working Environment:
 - Temperature -0C-40C
 - Relative Humidity 40% - 90%
- Display: 3 digits LED display with red words
- Dimension: 113 x 55 x 22mm (LxWxH), <80g
- Battery expectancy time: >200 times

1. Structure and Function Guide

Structure Guide

| | |
|---|--------------------|
| 1 | Mouthpiece |
| 2 | Vent hole |
| 3 | Display, LEDscreen |
| 4 | Start Button |
| 5 | Counter Button |
| 6 | Mini 5p USB Port |

Function Guide

Ready indicator light:

Low Battery display light:

Warning indicator light:

Ready for testing

Low battery indication

Excessive drinking

| | |
|---------------|----------------------|
| 0. 58 | Alcohol content |
| 20s Countdown | Warm-up Countdown |
| 10s Countdown | ready for blow |

2.Operation Instruction

- 2.1. Slide down the battery cover on the back side, insert 2 pcs AAA 1.5V batteries (not included) according to battery polarity indication and close it. Or use DC plug to directly connect with vehicle to use.
- 2.2. Press the start button last fort sec, the tester would be on with a buzzer brief ringing, Warm-up of 20 sec begin, the fig on the screen count down from 20-0, please wait now. After the countdown, tester will send out a ring, and the ready light will be on. Now the test could begin.
- 2.3. Have a deep breath before testing, and blow directly to gas entry until a buzzer ringing, "99" twinkles for 4 seconds on the screen.
- 2.4. Tester will send out a ring after testing process finished, the LED screen would give the readout of temporal alcohol concentration (If the testing result is lower than preset level, the warning light will be off; If the testing result is higher than preset level, the buzzer ringing continuously, at the same time the warning light will be red), this result would be keep on the screen for 10 sec, then cut off the power supply automatically with two buzzer ringing.
- 2.5. If a new testing needed, press the Start button again and repeat the above step 2.2 – 2.4.
- 2.6. Press the model button to choose the unit of the result you want.

3.Notification

- 3.1. Avoid any fall or strong shock.
- 3.2. If foreign gas of high concentration existed, may the tester won't work normally.

3.3. If testing is done under low voltage, certain error will be existed between the real value

3.4. To ensure the testing result, please wait 10 mins to take the testing after your drinking.

3.5. Please operate it strictly according to this instruction. If longtime not use, when you first begin anew testing, the first two times testing result may be unreliable, just begin testing from the third time, otherwise the result will also not be accurate.

3.6. Do not keep the tester in the environment of Corrosive gas (Chlorine etc) for using or storing, nor in other bad surroundings.

3.7. Under normal detections, sensor life could be more than two years.

3.8. After Long time use, there may be dirt on the tester, please use clean cloth to wipe off. Do not use any hard object or the solvent with any Corrosive ingredient.

Common Fault and Solving solutions

| Fault | Possible reason | Solution |
|------------------------------------|--------------------------------|--|
| NO display on LED | Incorrect battery installation | Insert battery correctly according to polarity |
| | Low battery seriously | New battery replacement |
| | Circuit fault | Contact distributor |
| No response to detection gas | not complete Warm-up | Waiting complete Warm-up |
| | Circuit fault | Contact distributor |
| FFF display then power off | Sensor fault | Contact distributor |
| Low voltage display then power off | Low battery seriously | New battery replacement |

The pictures are for reference only. The real design the products may differ from the pictures presented.

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